



Theophylline control (DAG-WT2707)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Concentration	$\geq 300 \mu\text{mol/L}$ (54.1 $\mu\text{g/mL}$)
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Keywords	Theophylline; 1,3-dimethylxanthine
----------	------------------------------------