



Isosorbide dinitrate control (DAG-WT2717)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	Unconjugated
Applications	Immunoassays
Format	Liquid or Lyophilized
Size	0.5 mL, 1.0 mL
Buffer	Synthetic serum
Storage	Store at -20°C

BACKGROUND

Introduction	Isosorbide dinitrate is a medication used for heart failure, esophageal spasms, and to treat and prevent chest pain from not enough blood flow to the heart. It has been found to be particularly useful in heart failure due to systolic dysfunction together with hydralazine. It is taken by mouth or under the tongue.
Keywords	Isosorbide dinitrate; ISDN; (3R,3aS,6S,6aS)-hexahydrofuro[3,2-b]furan-3,6-diyl dinitrate