



Recombinant *Felis domesticus* Fel d 4 (DAG-WT274)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	32 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	20mM Tris, 150mM NaCl, pH8.0
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle

BACKGROUND

Introduction	Fel d 4 is a lipocalin allergen produced by the cat, which binds IgE at relatively high frequency in cat-sensitive individuals. The allergen provides not only a means for investigating differences in the immune response to lipocalin allergens from that found for other mammalian species but also an important reagent for the diagnosis of cat allergy.
Keywords	Felis domesticus allergen; Fel d 4; Fel d 4.0101