



1-Aminohydantoin (AHD) [KLH] (DAG-WT2743)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Conjugate	KLH
-----------	-----

Applications	Immunogen
--------------	-----------

Format	Liquid
--------	--------

Size	1 mg
------	------

Buffer	PBS buffer
--------	------------

Storage	Store at -20°C
---------	----------------

BACKGROUND

Introduction	Aminohydantoin is a heterocyclic compound that belongs to the family of tetrazoles. It is an important intermediate in the synthesis of a variety of drugs and is commonly used as a precursor for the preparation of thiazoles like sulfonamides. Aminohydantoin is also used as a pH indicator, described by its characteristic color change from clear red to green or blue when dissolved in water.
--------------	---

Keywords	1-Aminoimidazolidine-2,4-dione; 1-aminohydantoin
----------	--