



Recombinant Human ICAM-1 (CD54) Protein [His] (DAG-WT2746)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 98% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Molecular Weight	66.2 kDa
Format	Lyophilized
Size	50 µg, 1 mg
Buffer	Lyophilized from sterile PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	ICAM-1 is a cell surface glycoprotein and an adhesion receptor that is best known for regulating leukocyte recruitment from circulation to sites of inflammation.
Keywords	Human intercellular adhesion molecule-1; ICAM-1/CD54; ICAM-1