



Recombinant *Aspergillus fumigatus* Asp f 1 (DAG-WT275)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	18 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	20mM Tris-HCl, 500mM NaCl, 8M urea, pH8.0
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle

BACKGROUND

Introduction	Asp f I is a major allergen produced by the mycelia of <i>Aspergillus fumigatus</i> . It is not present in spores and can be used as a specific marker for the detection of germination of this fungus.
Keywords	<i>Aspergillus fumigatus</i> allergen 1; Asp f 1; Asp f 1.0101