



Recombinant Human ADAMTS13 (aa 682 - 859) [His] (DAG-WT2752)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA, LF, WB
Molecular Weight	20.9 kDa
Format	Liquid or Lyophilized powder
Size	100 µg, 1 mg
Buffer	PBS, pH 7.4
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Cleaves the vWF multimers in plasma into smaller forms thereby controlling vWF-mediated platelet thrombus formation.
Keywords	ADAMTS13; A disintegrin and metalloproteinase with thrombospondin motifs 13; ADAM-TS 13; von Willebrand factor-cleaving protease

GENE INFORMATION

