



Recombinant Human APOE [His] (DAG-WT2754)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	35.3 kDa
Format	Lyophilized
Size	20 µg, 100 µg, 500 µg, 1 mg
Buffer	Lyophilized from 0.22 µm filtered solution in PBS (pH 7.4). Normally 5% trehalose is added as protectant before lyophilization
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Apolipoprotein E (APOE) is a class of apolipoprotein found in the chylomicron and IDLs that binds to a specific receptor on liver cells and peripheral cells. It is essential for the normal catabolism of triglyceride-rich lipoprotein constituents.
Keywords	Apolipoprotein E; APOE

GENE INFORMATION

UniProt ID

P02649
