



# Recombinant Human Beta-2-glycoprotein 1 Domain 1 (DAG-WT2778)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 85% as determined by SDS-PAGE
Conjugate	TBD
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	Tris-HCl buffer
Preservative	None
Storage	Store at -20°C

## BACKGROUND

Introduction	Binds to various kinds of negatively charged substances such as heparin, phospholipids, and dextran sulfate. May prevent activation of the intrinsic blood coagulation cascade by binding to phospholipids on the surface of damaged cells.
Keywords	ApoH; Beta-2-glycoprotein 1