



Recombinant Human RA33 protein (DAG-WT2780)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	Tris-HCl buffer
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Heterogeneous nuclear ribonucleoproteins A2/B1 is a protein that in humans is encoded by the HNRNPA2B1 gene.
Keywords	Heterogeneous nuclear ribonucleoprotein A2/B1; RA33; HNRNPA2B1