



Recombinant Human SSB.La (DAG-WT2783)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	Tris-HCl buffer
Preservative	None
Storage	Store at -20°C

BACKGROUND

Introduction	Binds to the 3' poly(U) terminus of nascent RNA polymerase III transcripts, protecting them from exonuclease digestion and facilitating their folding and maturation.
Keywords	SS-B/La; Lupus La protein; Sjogren syndrome type B antigen