



# Recombinant Human neuronal nuclear antigen (Ri) (DAG-WT2784)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 85% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	Tris-HCl buffer
Preservative	None
Storage	Store at -20°C

## BACKGROUND

Keywords	HuNNA; Neuronal nuclear antigen; Ri
----------	-------------------------------------