



# Recombinant Mus musculus Mus m 1 (DAG-WT279)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Purity</b>	> 90%
<b>Conjugate</b>	TBD
<b>Applications</b>	Immunoassays
<b>Molecular Weight</b>	17 kDa
<b>Format</b>	Liquid
<b>Concentration</b>	Batch dependent - please inquire should you have specific requirements
<b>Size</b>	1 mg
<b>Buffer</b>	20mM Tris, 150mM NaCl, pH8.0
<b>Preservative</b>	None
<b>Storage</b>	Store at -80°C. Avoid freeze / thaw cycle

## BACKGROUND

<b>Keywords</b>	Mus musculus allergen 1; Mus m 1; Lipocalin and urinary prealbumin; Mus m 1.0101; Mus m 1.0102
-----------------	--