



Recombinant *Ambrosia artemisiifolia* Amb a 2 (DAG-WT281)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	28 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	20mM Tris, 150mM NaCl, pH8.0
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle

BACKGROUND

Keywords	<i>Ambrosia artemisiifolia</i> allergen 2; Amb a 2
----------	--