



Recombinant *Chenopodium album* Che a 1 (DAG-WT284)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	35 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	50mM PB, pH7.8
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle

BACKGROUND

Introduction	Che a 1, a relevant allergen from chenopod pollen, is structurally related to the Ole e 1-like protein family, but exhibits significant differences on its polypeptide sequence that could explain its different antigenic behavior and limited cross-reactivity.
Keywords	<i>Chenopodium album</i> allergen 1; Che a 1; Ole e 1 homologue; Che a 1.0101