



# Recombinant *Chenopodium album* Che a 2 (DAG-WT285)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	15 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	20mM Tris-HCl, pH8.0
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle

## BACKGROUND

Keywords	<i>Chenopodium album</i> allergen 2; Che a 2; Profilin; Che a 2.0101
----------	--