



# Recombinant Humulus japonicus Hum j 1 (DAG-WT287)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	10 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	20mM Tris, 150mM NaCl, pH8.0
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle

## BACKGROUND

Keywords	Humulus japonicus allergen 1; Hum j 1; Hum j 1.0101
----------	---