



Recombinant Glycine max Gly m 5 (DAG-WT290)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	32 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle

BACKGROUND

Keywords	Glycine max allergen 5; Gly m 5; Beta-conglycinin
----------	---