



Recombinant Bos domesticus Bos d 2 (DAG-WT292)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	20 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	10mM PBS, pH7.4
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle

BACKGROUND

Keywords	Bos domesticus allergen 2; Bos d 2; Lipocalin
----------	---