



Recombinant Mangifera indica Man i 3 (DAG-WT296)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Purity	> 90%
Conjugate	TBD
Applications	Immunoassays
Molecular Weight	16 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	20mM Tris-HCl, pH8.0
Preservative	None
Storage	Store at -80°C. Avoid freeze / thaw cycle

BACKGROUND

Keywords	Mangifera indica allergen 3; Man i 3
----------	--------------------------------------