



Recombinant Zika Virus Envelope protein (a.a 40-400) [His] (DAG-WT302)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	Recombinant Zika envelope protein (strain MR766) is produced by E. coli and the target gene encoding a.a 40-400 is expressed with a His tag at the C-terminus.
Purity	> 95% , as determined by SDS-PAGE
Conjugate	His
Applications	N/A
Molecular Weight	38 kDa
Format	Liquid
Concentration	Batch dependent - please inquire should you have specific requirements
Size	1 mg
Buffer	PBS, with 25mM arginine
Preservative	None
Storage	2-8°C short term, -20°C long term

BACKGROUND

Introduction	Zika virus is an emerging disease that is spread by Aedes mosquitoes. The virus was first isolated in Central Africa, and has since been spread to South Asia and recently to South America. Zika virus can cause mild fever, rash, myalgia, arthralgia and headaches, with one in four infected individuals being asymptomatic. Due to similar symptoms Zika virus infected
---------------------	--

individuals can easily be mis-diagnosed as a dengue infection and vice-versa. In addition, Zika virus has been implicated in causing microcephaly through transmission in utero.

Keywords	Zika virus; ZIKV; Envelope protein
-----------------	------------------------------------
