



Recombinant ZIKV NS5 Protein [His] (DAG-WT310)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Antigen Description	A DNA sequence encoding the Zika virus (ZIKV) (strain ZikaSPH2015) NS5 (ALU33341.1) (Lys2772-Leu3423) was expressed with a polyhistidine tag at the C-terminus.
Purity	> 95% , as determined by SDS-PAGE
Conjugate	His
Applications	Immunoassays
Molecular Weight	76.7 kDa
Format	Lyophilized Powder
Size	1 mg
Buffer	20 mM PB,100 mM NaCl,50% glycerol, pH 7.0, 0.5 mM TCEP
Preservative	None
Storage	Store at -20°C to -80°C. Avoid multiple freeze/thaw cycles

BACKGROUND

Introduction	Zika virus is an emerging disease that is spread by Aedes mosquitoes. The virus was first isolated in Central Africa, and has since been spread to South Asia and recently to South America. Zika virus can cause mild fever, rash, myalgia, arthralgia and headaches, with one in four infected individuals being asymptomatic. Due to similar symptoms Zika virus infected individuals can easily be mis-diagnosed as a dengue infection and vice-versa. In addition, Zika virus has been implicated in causing microcephaly through transmission in utero.
---------------------	---

Keywords

Zika virus; ZIKV-NS5; Nonstructural protein

GENE INFORMATION

UniProt ID

ALU33341.1
