



# Native Cat Haptoglobin (DAG-WT3213)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Feline
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	ELISA
Format	Liquid
Size	100 µg, 1 mg
Buffer	10 mM sodium phosphate, pH 7.2, 150 mM NaCl
Preservative	0.1% sodium azide
Storage	Store at -20°C to -80°C

## BACKGROUND

Introduction	<p>Haptoglobin is a protein made by your liver. It attaches to a certain type of hemoglobin. Hemoglobin is a protein in your red blood cells that carries oxygen from your lungs to the rest of your body. Most hemoglobin is located inside red blood cells, but small amounts circulate in the bloodstream.</p>
Keywords	Haptoglobin; HP