



# Native thrombin-antithrombin complex (DAG-WT3235)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Calibrator, Standard
Format	Liquid
Size	100 µg
Buffer	10mM PBS, 20mM Tris, Protein stabilizer
Preservative	0.1% Proclin300
Storage	Store at -20°C

## BACKGROUND

Introduction	Thrombin–antithrombin complex (TAT) is a protein complex of thrombin and antithrombin. It is a marker of net activation of coagulation.
Keywords	Thrombin–antithrombin complex; TAT