



# Recombinant Human HBP (DAG-WT3236)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 90% as determined by SDS-PAGE
Conjugate	His
Applications	ELISA
Format	Liquid
Size	1 mg
Buffer	Tris/PBS-based buffer
Storage	Store at -20°C to -80°C

## BACKGROUND

Introduction	Heparin-binding protein (HBP, also known as azurocidin or CAP37), stored by neutrophils, is released early upon neutrophil adhesion and extravasation. Bacterial products induce HBP release and HBP is the key mediator of neutrophil-dependent vascular permeability increases. Elevated plasma levels of HBP are a promising early biomarker of severe sepsis in the emergency department and correlate with the development of septic shock.
Keywords	Heparin-binding protein; HBP