



# Recombinant Human MPO (DAG-WT3245)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

|              |                                 |
|--------------|---------------------------------|
| Species      | Human                           |
| Purity       | > 90% as determined by SDS-PAGE |
| Conjugate    | His                             |
| Applications | ELISA                           |
| Format       | Liquid                          |
| Size         | 1 mg                            |
| Buffer       | Tris/PBS-based buffer           |
| Storage      | Store at -20°C to -80°C         |

## BACKGROUND

|          |                      |
|----------|----------------------|
| Keywords | Myeloperoxidase; MPO |
|----------|----------------------|