



## Recombinant Human B2-MG (DAG-WT3252)

This product is for research use only and is not intended for diagnostic use.

### PRODUCT INFORMATION

|              |                                 |
|--------------|---------------------------------|
| Species      | Human                           |
| Purity       | > 90% as determined by SDS-PAGE |
| Conjugate    | His                             |
| Applications | ELISA                           |
| Format       | Liquid                          |
| Size         | 1 mg                            |
| Buffer       | Tris/PBS-based buffer           |
| Storage      | Store at -20°C to -80°C         |

### BACKGROUND

|          |                                  |
|----------|----------------------------------|
| Keywords | B2-MG; Beta 2-Microglobulin; B2M |
|----------|----------------------------------|