



# Recombinant Human Coagulation Factor VIII (DAG-WT3259)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 97% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	Immunoassays
Reconstitution	Reconstitute 250IU lyophilized Factor VIII in 5ml sterile water
Format	Lyophilized
Size	250 IU, 500 IU, 1000 IU
Buffer	Lyophilized from a solution containing 8mg Tween-80, 112mM NaCl, 40mg Mannitol, 10mg Trehalose, 1ng VWF and 4.2mM CaCl <sub>2</sub>
Storage	Store at -20°C

## BACKGROUND

Keywords	Factor VIII; FVIII; Anti-hemophilic factor (AHF)
----------	--