



Recombinant Human Coagulation Factor VIIa (DAG-WT3261)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 98% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	Immunoassays
Molecular Weight	50 kDa
Format	Lyophilized
Size	10 µg, 50 µg, 1 mg
Buffer	Lyophilized after from a sterile solution containing 10mg sucrose pH6
Storage	Store at -20°C

BACKGROUND

Keywords	Coagulation Factor VIIa; FVIIa; rFVIIa
----------	--