



Recombinant Human G-CSF (aa30-204) (DAG-WT3264)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	Unconjugated
Applications	Immunoassays
Molecular Weight	17.8 kDa
Format	Liquid
Size	2 µg, 10 µg, 50 µg, 100 µg
Buffer	PBS, pH 7.4
Storage	Store at -20°C

BACKGROUND

Keywords CSF3; Granulocyte colony-stimulating factor; G-CSF; Pluripoietin; Filgrastim; Lenograstim

GENE INFORMATION

UniProt ID P09919-2