



Recombinant Human G-CSF [His, MBP] (DAG-WT3265)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	His, MBP
Applications	ELISA, WB
Molecular Weight	61.5 kDa
Format	Liquid
Size	20 µg
Buffer	0.15 M Phosphate buffered saline, pH 7.4
Preservative	None
Storage	Store at -20°C to -80°C

BACKGROUND

Keywords CSF3; Granulocyte colony-stimulating factor; G-CSF; Pluripoietin; Filgrastim; Lenograstim

GENE INFORMATION

UniProt ID P09919