



Native Human Factor VII (DAG-WT3274)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	50 kDa
Format	Liquid or Lyophilized powder
Size	100 µg, 1 mg
Buffer	20 mM Tris-HCl, 100mM NaCl, 1 mM Benzamidine, pH 7.4
Storage	Store at -80°C

BACKGROUND

Keywords	Factor VII; Coagulation factor VII; EC 3.4.21.21; Proconvertin
----------	--