



Native Human Activated Protein C (DAG-WT3276)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

| | |
|------------------|--|
| Product Overview | The Activated Human Protein C was activated from homogeneous protein C by activation with purified alaphthrombin. The alpha-thrombin is removed after activation |
| Species | Human |
| Conjugate | N/A |
| Applications | Immunoassays |
| Molecular Weight | 56 kDa |
| Format | Liquid |
| Size | 100 µg, 1 mg |
| Buffer | 20 mM Tris-HCl, 100 mM NaCl, pH 7.4 |
| Storage | Store at -80°C |

BACKGROUND

| | |
|----------|--------------------------|
| Keywords | Activated Protein C; APC |
|----------|--------------------------|