



Native Bovine Activated Protein C (DAG-WT3278)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The Activated Bovine Protein C was activated from homogeneous protein C by activation with purified alaphthrombin
Species	Bovine
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	52.7 kDa
Format	Liquid
Size	100 µg
Buffer	20 mM Tris-HCl, 100 mM NaCl, pH 7.4
Storage	Store at -80°C

BACKGROUND

Keywords	Activated Protein C; APC
-----------------	--------------------------