



Native Human Factor IXa Beta (DAG-WT3280)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	45 kDa
Format	Liquid
Size	500 µg
Buffer	50mM MES, 0.6M NaCl, pH 6.0
Storage	Store at -80°C

BACKGROUND

Keywords	Factor IX; Christmas factor; EC 3.4.21.2
----------	--