



Native Bovine Factor IX (DAG-WT3281)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Bovine
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	55.4 kDa
Format	Liquid
Size	500 µg
Buffer	20 mM Tris-HCl, 100mM NaCl, 1 mM Benzamidine, pH 7.4
Storage	Store at -80°C

BACKGROUND

Keywords	Factor IX; Christmas factor; EC 3.4.21.2
----------	--