



Native Bovine Factor X (DAG-WT3286)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Bovine
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	55 kDa
Format	Liquid
Size	1 mg
Buffer	20 mM Tris-HCl, 100mM NaCl, 1 mM Benzamidine, pH 7.4
Storage	Store at -80°C

BACKGROUND

Keywords	Factor X; F10; FX; FXA; Coagulation factor X
----------	--