



Native Bovine Factor Xa (DAG-WT3287)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Bovine
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	45.3 kDa
Format	Liquid
Size	1 mg
Buffer	20 mM Tris-HCl, 0.7 M NaCl, pH 7.4
Storage	Store at -80°C

BACKGROUND

Keywords	Factor Xa; Coagulation factor X; Thrombokinase; Stuart-Prower factor
----------	--