



Native Human Factor XIa (DAG-WT3294)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	160 kDa
Format	Liquid
Size	100 µg, 1 mg
Buffer	4 mM Sodium Acetate-HCl, 0.15 M NaCl, pH 5.3
Storage	Store at -80°C

BACKGROUND

Keywords	Factor XI; Plasma thromboplastin antecedent; F11
----------	--------------------------------------------------