



# Native Human Factor XIIa beta (DAG-WT3298)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	30 kDa
Format	Liquid
Size	100 µg, 1 mg
Buffer	4 mM Sodium Acetate-HCl, 0.15 M NaCl, pH 5.3
Storage	Store at -80°C

## BACKGROUND

Keywords	Factor XII; Coagulation factor XII; Hageman factor
----------	--