



Native Human Factor XIII (DAG-WT3299)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	320 kDa
Format	Liquid
Size	250 µg, 1 mg
Buffer	50 mM Tris-HCl, 0.1 M NaCl, 1 mM EDTA, 10u/mL aprotinin, 20% glycerol, pH 7.5
Storage	Store at -80°C

BACKGROUND

Keywords	Factor XIII; Factor IXa; Fibrin stabilizing factor
----------	--