



Native Human Factor XIIIa (DAG-WT3300)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	312 kDa
Format	Liquid
Size	250 µg, 1 mg
Buffer	0.05 mM EDT, 50% glycerol
Storage	Store at -80°C

BACKGROUND

Keywords	Factor XIII; Factor IXa; Fibrin stabilizing factor
----------	--