



Native Murine Factor XIII (DAG-WT3301)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Murine
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Format	Liquid
Size	100 µg
Buffer	50 mM Tris, 100 mM NaCl, 1 mM EDTA, 10 u/mL aprotinin, 20% glycerol, pH 7.5
Storage	Store at -80°C

BACKGROUND

Keywords	Factor XIII; Factor IXa; Fibrin stabilizing factor
----------	--