



Native Human Factor VIIa (DAG-WT3302)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The Human Factor XIIa was activated from homogeneous Human Factor XII using Dextran Sulfate. This Dextran Sulfate was removed after activation.
Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	50 kDa
Format	Liquid
Size	100 µg, 1 mg
Buffer	20 mM Tris-HCl, 0.1 M NaCl, pH 7.4
Storage	Store at -80°C

BACKGROUND

Keywords	Coagulation factor Vii; FVII; rFVII
-----------------	-------------------------------------