



Native Platelet Factor 4 (DAG-WT3303)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	29 kDa
Format	Liquid
Size	100 µg
Buffer	140 mM Tris-HCl, 1.05 M NaCl, pH 7.4
Preservative	0.02% NaN ₃
Storage	Store at 2 to 8°C

BACKGROUND

Keywords	Factor IV; Calcium
----------	--------------------