



Native High Molecular Weight Kininogen (Single Chain) (DAG-WT3306)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	120 kDa
Format	Liquid
Size	1 mg
Buffer	4 mM Sodium Acetate-HCl, 0.15 M NaCl, pH 5.3
Storage	Store at -80°C

BACKGROUND

Keywords	Kininogen; KNG
----------	----------------