



Native Bovine Plasminogen (DAG-WT3309)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Bovine
Purity	> 90% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	90 kDa
Format	Liquid
Size	1 mg
Buffer	50 mM Tris-HCl, 0.1 M NaCl, pH 7.4
Storage	Store at -80°C

BACKGROUND

Keywords	PLG; Plasminogen; Plasmin; HAE4
----------	---------------------------------