



Native Human γ Thrombin (DAG-WT3311)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Species	Human
Purity	> 95% as determined by SDS-PAGE
Conjugate	N/A
Applications	Immunoassays
Molecular Weight	12 kDa
Format	Liquid
Size	1 mg
Buffer	10 mM Sodium Acetate, 0.75 M NaCl, pH 6.0
Storage	Store at -80°C

BACKGROUND

Keywords	Thrombin; Gamma thrombin; Coagulation factor II
----------	-------------------------------------------------