



Native Canine Thrombin (DAG-WT3312)

This product is for research use only and is not intended for diagnostic use.

PRODUCT INFORMATION

Product Overview	The Canine Thrombin was activated from canine prothrombin by activation with venom. This venom is removed after activation
Species	Canine
Conjugate	N/A
Applications	Immunoassays
Format	Liquid
Size	100 µg
Buffer	50 mM sodium citrate, 0.2 M NaCl, pH 6.5
Preservative	0.1% PEG-8000
Storage	Store at -80°C

BACKGROUND

Keywords	Thrombin; alpha thrombin; Coagulation factor II
-----------------	---