



# Native Murine Thrombin (DAG-WT3313)

This product is for research use only and is not intended for diagnostic use.

## PRODUCT INFORMATION

<b>Product Overview</b>	The Murine Thrombin was activated from homogeneous Murine Prothrombin by activation with Saw Scaled Viper Venom. The venom was removed after activation
<b>Species</b>	Murine
<b>Conjugate</b>	N/A
<b>Applications</b>	Immunoassays
<b>Format</b>	Liquid
<b>Size</b>	50 µg
<b>Buffer</b>	50 mM sodium citrate, 0.2 M NaCl, pH 6.5
<b>Preservative</b>	0.1% PEG-8000
<b>Storage</b>	Store at -80°C

## BACKGROUND

<b>Keywords</b>	Thrombin; alpha thrombin; Coagulation factor II
-----------------	---